
KL 2500 LED User Manual

DE | EN | FR | ES | IT

SCHOTT
glass made of ideas

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Instrument overview

❶ Electronic light intensity setting

❷ On/off switch

❸ Locking screw

❹ Light guide socket

❺ Carrying handle

❻ 24V connection

❼ Ventilation grid

❽ Filter slide

❾ Air vents

❿ Information label

⓫ Indicator light

⓬ LCD display

⓭ USB connection socket

⓮ Foot switch connection socket

Content

1. Important information	14
2. Operation	16
2.1 Light guide connection	16
2.2 Power connection	16
2.3 Start-up procedure	17
2.4 Light intensity setting	17
2.5 Filter slide	18
2.6 Control via USB	19
2.7 Control via foot switch	19
3. Maintenance	20
4. Troubleshooting	20
5. Technical data	21

1. Important information

Symbols used

Symbol Meaning

Warning of danger
(caution, obey documentation)

LED-radiation (caution, do not stare into beam!)

Intended use

The KL 2500 LED cold light source is intended for use in industrial and laboratory settings¹. Cold light sources are used to illuminate objects in stereomicroscopy. The visible and intense light from the cold light source is guided to the object via light guides. Only SCHOTT light guides and SCHOTT accessories ensure proper function, safety and optimum light output.

Safety information

Please read and follow these instructions carefully. Ensure that every operator of the device has quick access to these operating instructions.

The instrument's safety cannot be guaranteed if the instructions are not followed. The manufacturer is not liable for any damage caused by non-observance of these operating instructions.

- Operate the light source only in dry conditions and not in hazardous areas.
- When connecting the device, ensure a safe connection.
- Operate the cold light source only with the voltage specified on the information label.
- Never look directly into the light emission from the light source or the light guide during operation (risk of eye injuries)! The KL 2500 LED emits visible light of high intensity (according to the EN 62471:2008 standard, risk class 2).
- Never cover the open light guide receptacle or the light guide output (fire hazard)!
- Never cover the open light guide receptacle or the light guide outlet with your hand or other parts of your body (risk of burns)!

¹ You can obtain an individual spectrum of your light source at info.microscopy@schott.com

- When illuminating heat-sensitive or flammable light-absorbing objects, pay particular attention to selecting the light guide spacing and intensity so that no thermal damage occurs to the object. Since light-absorbing materials have the physical property of converting incident light into heat, damage can occur to heat-sensitive or flammable light-absorbing materials.
- All light guide outputs must always be at a safe distance – at least 10 cm – from heat-sensitive or flammable light-absorbing materials when the light source is switched on (avoid possible fire hazard).
- Reduce the brightness and duration of illumination of biological tissue to avoid unnecessary exposure of biological tissue to visible light illumination.
- Keep all ventilation openings free. In case of insufficient cooling, an integrated thermal switch will regulate the light intensity of the light source.
- Safe disconnection from the power supply is only achieved by pulling out the plug.
- Do not open or disassemble the device. Refrain from making technical modifications to the device.
- Allow the filter slider and insert filter to cool down before removing the insert filter. The filter slider must always be in a detent position when operating the light source (see "Filter slide").
- Allow the light guide to cool down before replacing it. We recommend that you switch off the light source for about five minutes before replacing the light guide. Do not touch the light guide outlet with your hand or other parts of your body (risk of burns)!

2. Operation

2.1 Light guide connection

First, open the locking screw of the light guide holder by turning it counter-clockwise. Then push the light guide in as far as it will go and tighten the locking screw again by turning it clockwise.

When inserting light guides with a guide pin, make sure that it aligns in one of the two guide slots.

2.2 Power connection

Insert the hollow plug of the enclosed power supply unit into the socket on the light source.

Ensure that the KL 2500 LED cold light source is operated only in the specified voltage range and with the enclosed power supply unit.

Plug the three-pole power cable of type H05VV-F3G0.75 (available as an accessory) into the connection socket of the enclosed power supply unit provided for this purpose. When replacing the power cord, use only power cords of the same type.

You can switch the KL 2500 LED on and off by pressing the main switch. If the device is switched on, the green control lamp and the LCD display light up. To disconnect the power supply, pull out the plug!

2.3 Start-up procedure

Adjust the light intensity by turning the light intensity control. The light intensity is shown in the LCD display in values from 0% to 100%.

2.4 Light intensity setting

You can adjust the light intensity in three sensitivity levels:

1. Standard: in 1% steps
2. Coarse: in 5% steps. "Coarse" is shown in the display.
3. Fine: in 0.1% steps. "Fine" is shown in the display.

Briefly press the light intensity control to switch between the sensitivity levels.

Attention: The device is not switched off in the 0.0% control position, therefore the green control lamp lights up. To switch off the KL 2500 LED, press the main switch.

The percentage display in the LCD display corresponds to the respective current consumption of the LEDs. The control is linear. However, the radiant power of the LEDs is not proportional to the current consumption. The intensity change in the control range at more than 50.0% is therefore lower than in the lower control range.

Press the light intensity control for more than two seconds to manually save the displayed brightness. "Saved" then appears in the LCD display. If a light intensity has been saved, the KL 2500 LED light source lights up with this saved light intensity in the "Standard (1%)" sensitivity level after it has been switched off and on again.

2.5 Filter slide


The KL 2500 LED has a filter slide that can be equipped with an insert filter. Various filters are available as accessories.

When inserting the filter slide into the cold light source, there are two detent positions which are signaled by a clicking sound. The filter slide must be in one of the two end positions or the detent position when operating the light source. If the filter slide is in an intermediate position, it may damage the light source.


Insert filter

Make sure that the filter slide has cooled down before inserting a filter. Wear lint-free gloves to load the filter slide.


1. Pull out the filter slide completely and turn it so that the labeling is legible.

2. Insert the insert filter from the top right behind the two retaining lugs marked in black.


3. Click the insert filter firmly into place by applying slight pressure at the top right. The filter now lies flat in the holder.

To remove, bend the retaining clip slightly to the left and remove the filter.

Remove filter


If you want to operate the light source without a filter for a short time, pull the slider out only to the first detent point. In this position, the filter is still in the light source, but no longer in the beam path.

Your KL 2500 LED has a USB interface and can be controlled via a computer, by using the enclosed USB connection cable. All information such as port settings, formats, control commands and error codes required for integrating the KL 2500 LED into individual control software is described in a separate information sheet.

2.6 Control via USB

In addition, SCHOTT offers software for demonstration purposes. The overview of the protocol commands and the demo software can be downloaded from the SCHOTT website <https://schott.com/products/kl-fiber-optic-light-sources> download.

First, install the appropriate driver so that the KL 2500 LED can communicate with your computer. To do this, connect the light source to your computer using the enclosed USB connection cable and switch on the light source. The operating system will prompt you to install the new USB device. Follow the prompts of your operating system.

The illumination of the KL 2500 LED can be switched on and off using an optionally available foot switch. To do this, connect the foot switch to the corresponding connection socket on the back of the light source.

2.7 Control via foot switch

If the lighting is switched off by means of the foot switch, this is indicated on the display "SHUTTER". In addition, the light intensity control is locked. Press the foot switch again to unlock the light intensity control and set the original brightness.

3. Maintenance

Your KL 2500 LED is maintenance-free.

Use a soft dry cloth or commercially available plastic cleaning cloths for external cleaning of the device.

4. Troubleshooting

If you are unable to rectify the problem using the measures described below, please contact your specialist dealer or the nearest SCHOTT agency.


Type of fault	Possible cause	Remedy
1. Lamp out	<ul style="list-style-type: none"> • Instrument not switched on • Plug not in socket • No mains electricity voltage 	<ul style="list-style-type: none"> • Switch instrument on • Plug in socket • Check mains voltage
2. Light intensity dropping	<ul style="list-style-type: none"> • Electronics overheated 	<ul style="list-style-type: none"> • Switch off instrument. Ensure adequate cooling, switch on again after prolonged cooling time.
3. Lamp off and display shows „SHUTTER“	<ul style="list-style-type: none"> • Accidental electrostatic discharge of the connection socket of the foot switch 	<ul style="list-style-type: none"> • Switch the device off and on again using the main switch

5. Technical data

Properties	Values	
General information		
Type description		KL 2500 LED
Dimensions (W x D x H)	mm	approx. 114 x 231 x 137
Weight	kg	approx. 2,45
Cooling		Axial fan
Ambient temperature*	°C	+ 5 ... + 40
Relative air humidity*	%	up to 31°C ambient temperature: 85 31°C to 40°C Ambient temperature: decreasing linearly to 75%
Air pressure*	hPa	800 ... 1060
Transport and storage		
• Temperature	°C	-20 ... +70
• Rel. air humidity	%	10 ... 95 (non-condensing)
• Air pressure	hPa	500 ... 1200
Contamination level		2
IP-Protection class		IP20
Elektrical information		
Operating voltage, frequency		100 - 240 V ~ 50 - 60 Hz
Nominal voltage light source	V	24 DC
Rated current light source	A	2,1
Power consumption	VA	max. 51
Protection class power supply		I
Protection class light source		III
Oversvoltage category		II
Lamp type		High Power LED
Average LED lifetime		
Setting 100%	h	50.000 (luminous flux reduced to 70%)
Lighting information		
Active area of the light guide	mm	9
Total luminous flux at light guide exit (SCHOTT light guide Ø 9 mm, typical value)		
100% setting	lm	1,100
Color temperature	K	Approx. 5,600
Light exit angle ($2\alpha_{\text{eff}}$)	°	Approx. 70
Conformity		CE, UKCA (power supply CE, UKCA, UL, PSE)
EMC emission classe		B

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WEEE declaration

Your SCHOTT product was produced and developed with high quality materials and components. The symbol indicates that electrical and electronic devices must be separated from domestic waste and appropriately disposed of after useful life. SCHOTT AG Lighting and Imaging has arranged a waste management system for recycling. Please use this system for removal and help to protect the environment we live in.

Further information regarding our waste management system, please refer to www.schott.com/lightingimaging/recycle

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SCHOTT North America, Inc.