




MicroscopeWorld

ZEISS ZEN Lite Quick Start Guide

How to install ZEN Lite and view and save images.

Installation

First, connect your power cable and USB connection cable to your camera and your USB port on your computer/laptop.



Plug the other end of this cable to a USB 3 port on your computer.


Make sure to plug the other end of this cable to a power source. A blue light should come on at the top of the camera, indicating that it is on.

Next, download the software and drivers for your camera.

Go to <https://portal.zeiss.com/download-center/>

**You will need to register with Zeiss to be able to download the software.

After you've registered, click on the "Download Center" option on the left hand side of the page.



Installation

Once in the download center, click on the “ZEN/ ZEN light” option.



ZEN / ZEN lite

Latest version: 3.6


Digital Imaging Suite for light and confocal microscopy

[Show all files](#)

Click on the blue Download Installer link to begin the download process.

This ZEN product is installed through ZEISS Microscopy Installer

[Download Installer →](#)



ZEN 3.6 - Full Version

License needed

Machine Compatibility: Widefield Systems, LSM 800 for materials, LSM 900 for materials, LSM 800, LSM 900, LSM 980, ZEN Correlative Array Tomography (CAT), Celldiscoverer 7, Lightsheet 7, Lightsheet.Z1

OS Compatibility: Microsoft Windows 10 x64 Professional or Enterprise

Released: 05/19/2022

Release Highlights

- **Intellesis Denoising (AI Toolkit):** A deep-learning based method to learn “the noise” from the image itself (no labels or ground truth need) in order to remove the noise. User friendly training tools, annotation need not be provided, supports multi-dimensional data in ZEN (Channels, Z-stacks, Tiles).
- **ZEN Connect (Connect Toolkit):** ZEN Connect is supported by a manual 3D alignment with display of two volumes from all three angles and in 3D. Shifting the selected volume and changing opacity and LUT. Alignment with the help of a dedicated alignment tab.
- **Direct Processing:** Setup optimized for using only one workstation. New functions available: Linear Unmixing

[Read More](#)



Installation

Click on the blue Download button to install the ZEISS Microscopy Installer.

[Home](#) / [Download Center](#) / [ZEISS Microscopy Software](#) / [ZEISS Microscopy Installer](#)

ZEISS Microscopy Installer

[← Back to overview](#) Interested in upcoming releases? [🔔 Notify Me](#)



ZEISS Microscopy Installer

OS Compatibility: Microsoft Windows 10 x64 Professional or Enterprise
Released: 08/25/2022
Size: 93.5 MB

[Download](#)

Release Highlights

The ZEISS Microscopy Installer application allows you to install, uninstall, and update ZEISS Microscopy software products directly from the cloud.

For the new features of Version 1.5, read the [Release Notes](#).

For more information, see [User Manual](#).

It supports the following products (also including components, such as MTR2011, Python, Word Add-In):


[Read More](#)

Important Notes

- Make sure that you are entitled to use a selected software product before installing it through the ZEISS Microscopy Installer. ZEISS Microscopy Installer will allow you to install any software and software version from its product portfolio, even those that you are not licensed to use. These products will not function without a corresponding license.
- Do not upgrade your software product to a higher software version if you are unsure whether your license is valid for that version of the software. Installing an unlicensed version of the software might render your microscopy system or workstation nonfunctional.


[Read More](#)

Once the ZEISS Microscopy Installer installation is complete, the icon below will appear on your desktop. Double click the icon to open the installer.



Installation

The Microscopy Installer will pull up the screen below. Highlight the desired program to install (ZEN/ ZEN lite). Make sure that the “I agree to the Terms and Conditions” button in the lower left is selected, and then select “Log in to install”.



The ZEISS installer will have you fill in a form your information, which will generate an email so that ZEISS can grant you approval to download the software. After receiving the email, click on the emailed link to complete the installation.

The form is titled 'ZEISS ID' and 'Your access to ZEISS Microscopy Installer'. It shows a progress indicator with two steps: '1 Create your ZEISS ID account' and '2 Your access to ZEISS Microscopy Installer'. Below this, a message states: 'Thank you. Your ZEISS ID account has been created. For your access to ZEISS Microscopy Installer we need some more information from you.'

The form fields are:

- Salutation: -
- First Name*: Charisa
- Last Name*: Clarkson
- Email*: [empty]
- Language: [empty]

The ZEISS logo and a microscope icon are displayed on the right side of the form.




Image Viewing with ZEN lite

- Make sure that the camera and microscope are on prior to opening the ZEN software.
- Double click the ZEN icon on your desktop to open the program.



The following screen will open. By default, this welcome screen will always appear. If you prefer not to see this screen every time you open the program, check the box labeled “do not show this at startup”




You can click on the appropriate button if you would like to open an existing image, acquire a new image, or create a new project, or you can click the x at in the right hand corner of the box to close it.



Viewing a live image

- First, make sure that your camera and microscope are on. By default, the “Camera” tab should be selected. (Click on it if not.)
- Place a sample on the microscope and use the eyepieces to bring the sample into focus.



If the camera is connected properly, the “Active camera” should indicate the type of camera connected to your microscope. Otherwise, select the appropriate camera from the drop down menu. Note that the camera must be powered up prior to opening the software in order for ZEN to connect to it.

In the Microscope components, make sure to check “Show all” box. Change the microscope, objective, and other parameters to match your microscope set-up.
NOTE: The “Total Magnification” value will not match the true microscope magnification.


On the “camera” menu, checking the Auto-exposure box may help to bring the image into focus. Try selecting and deselecting this setting if you are having difficulties getting the camera to focus.

To perform a white balance, which helps ensure proper sample coloration, place a white piece of paper under the objective, and hit the “Auto” button under “White Balance.”



Capturing an image


- With your sample under the microscope objective, click the “Live” button. After clicking the Live button, the live button will turn red. Your screen should show an image of your sample that is continuously updating.
- The image may be adjusted by manipulating the settings in the Camera box.



Capturing an image



- Moving the tabs along the bottom on the graph on the Display tab can help to balance the colors of your image.
- Clicking "Reset" will restore the original settings if needed.




To take a picture of the sample, click "Snap".




Changing the default location for saving an image

- By default, the ZEN Program will save files on your computer's hard drive. To change or to view where the images are being saved, open the options tab under the Tools menu along the top left of the screen.
- Select the "Saving" path to view where images are being saved. To change this path, click the "..." button to bring up dialog box where you can select a different default saving location.
- You can check the "Auto Save after Snap" if you would like to save the image to the default location each time the "Snap" button is clicked.



Saving an image manually

- To save an image, click “Save” to save the image to the default location set up by the ZEN program. (Refer to the previous page for instructions to change this location.) By default, the image will be saved as a .czi file type.
- If another file type is desired (i.e. .jpeg), from the File menu in the upper left hand part of the screen, select “Save as with Options”. A dialog box will appear. Change the file format in the drop down menu of the dialog box.





MicroscopeWorld

www.MicroscopeWorld.com
6122 Innovation Way
Carlsbad, CA 92009
800-942-0528
info@microscopeworld.com

