

# CS Series Centrifuges

CS1385 / CS2385 / CS1465 / CS2465

Precision. Performance. Reliability.

CENTRIFUGES

The CS Series centrifuges deliver reliable, compact performance for clinical and laboratory applications. Featuring fixed-angle and horizontal rotor options, variable or fixed speed control, and a universal power supply, the CS1385, CS2385, CS1465, and CS2465 are engineered for precision sample processing in space-conscious environments.

ROTOR & CAPACITY				
Specification	CS1385	CS2385	CS1465	CS2465
Rotor Type	Fixed Angle		Horizontal	
Tubes	8		6	
Tube Positions	8×1.5/2 mL		6×15 mL	

SPEED & PERFORMANCE				
Specification	CS1385	CS2385	CS1465	CS2465
Speed	1,628 RPM (Fixed)		3,400 RPM (Fixed)	500–4,000 RPM (Variable)
Max RCF	1,628 g		1,712 g	

CONTROLS & DISPLAY				
Specification	CS1385	CS2385	CS1465	CS2465
Display	Segment LCD			
Operation	Keys			
Adjust Speed During Run	—			
Preset Groups	3 Groups			

ELECTRICAL & PHYSICAL				
Specification	CS1385	CS2385	CS1465	CS2465
Power Supply	AC 100–240V, 50/60Hz, 120W			
Dimensions (W×D×H)	12.4×12.4×15.5 in / 316×394×250 mm			
Weight	11 lb / 5 kg			

Specifications subject to change without notice.

© 2026 MPD Scientific. All rights reserved.

