

RDT3510

10-Position Fixed-Speed Single-Layer Mixer

Reliable Everyday Mixing, Right-Sized for Your Lab



10

POSITIONS

24

RPM

8W

POWER

1

LAYER

10 Positions

Single-Layer Design

24 RPM

Fixed End-Over-End

Type A Roller

Standard Length

8W / 110V 60Hz

Energy Efficient

WHERE THE RDT3510 PERFORMS BEST

Clinical Laboratories

Consistent end-over-end rotation for anticoagulated blood tubes, plasma samples, and reagents ahead of hematology and coagulation analyzers.

Incubators & Cold Rooms

Compact single-layer design and 8W power draw make it the ideal mixer for temperature-controlled environments where space and heat output matter.

Research Facilities

Gentle long-duration agitation for cell suspensions, buffer solutions, and biological specimens.

Diagnostic Laboratories

Everyday reagent homogenization and sample prep workflows where a simple, low-maintenance, consistent mixer is exactly what's needed.

FULL SPECIFICATIONS

Model	RDT3510	Roller Type	Type A (Standard Length)	Speed	24 RPM
Type	Fixed-Speed Mixer	Positions	10	Rotation	End-Over-End
Configuration	Single-Layer	UPC	00850084643323	Power	AC 110V, 60Hz, 8W