

FIXED-SPEED MIXER

Double the positions. Same footprint.

More samples, less space.



20
POSITIONS

24 RPM
FIXED SPEED

Type A
FORMAT



APPLICATIONS

The RDT3720 is built for labs that need to process high sample volumes without taking up additional bench space. Its 20-position dual-layer design is ideal for clinical and research environments handling blood, serum, reagents, or biological samples requiring continuous gentle mixing. The compact Type A footprint fits standard incubators and cold rooms, making it a practical choice for temperature-sensitive protocols across hematology, immunology, and general sample preparation workflows.

KEY FEATURES

- Double-layer design — 20 positions in the space of 10
- Fixed-speed end-over-end rotation at 24 RPM
- Type A compact footprint for incubators & cold rooms
- Uniform mixing across all positions simultaneously
- Ideal for blood, reagent & biological sample processing
- CE certified, 1-year warranty

TECHNICAL SPECIFICATIONS

Product Type	Tube Rotator
Speed	24 RPM (Fixed)
Positions	20 (2 layers x 10)
Format	Type A
Power Input	AC 100-240V / 0.6A
Power Output	DC 24V / 1A
Certification	CE
Warranty	1 Year
Origin	China

COMPATIBLE WITH 50+ TUBE FORMATS

Works with a wide range of standard laboratory tubes including EDTA, serum, plasma, and microcentrifuge tubes. See full compatibility list at mpdscientific.com.