

# RDT3720 Fixed-Speed Mixer

Fixed-speed double-layer mixer — 20 positions, Type A

SKU: RDT3720 • UPC: 00850084643347

● ACTIVE

TYPE A

The RDT3720 is a fixed-speed double-layer tube rotator that doubles the capacity of the RDT3510 by stacking two rows of 10 tube positions — totaling 20 positions — within the same compact Type A footprint. The dual-layer design maximizes throughput without requiring additional bench or incubator space, making it an efficient choice for labs processing high sample volumes. Consistent fixed-speed end-over-end rotation ensures reliable, gentle mixing performance across all 20 positions simultaneously.

Positions

**20**

2 layers x 10

Speed

**24 RPM**

Fixed — end-over-end

Layers

**2**

Double-layer design

KEY FEATURES	
<ul style="list-style-type: none"> <li>• Double-layer design with 20 positions — twice the capacity of single-layer models</li> </ul>	
<ul style="list-style-type: none"> <li>• Fixed-speed end-over-end rotation for uniform mixing across all positions</li> </ul>	
<ul style="list-style-type: none"> <li>• Type A format — compact footprint fits standard incubators and cold rooms</li> </ul>	
<ul style="list-style-type: none"> <li>• Maximizes throughput without additional bench space</li> </ul>	
<ul style="list-style-type: none"> <li>• Ideal for high-volume blood, reagent, and biological sample processing</li> </ul>	

TECHNICAL SPECIFICATIONS	
Product type	<b>Tube Rotator</b>
Speed type	<b>Fixed-Speed</b>
Speed	<b>24 RPM</b>
Total positions	<b>20</b>
Layers	<b>2</b>
Layer style	<b>Double-Layer</b>
Format	<b>Type A</b>
Origin	<b>China</b>

POWER & ELECTRICAL	
Included accessory	<b>Power Adapter</b>
Input	<b>AC 100-240V / 0.6A</b>
Output	<b>DC 24V / 1A</b>

COMPLIANCE & WARRANTY	
Certification	<b>CE</b>
Warranty	<b>1 year</b>

**COMPATIBLE PRODUCTS (50 items)**

BLAUNUFT0100ULU	BLAP11065U	BLAP11073U	BLAP11077U	BLAUNUFT0300UL-LRU	BLAUNUFT0100UL-LRU
BLAUNUFT0200UL-LRU	BLAUNUFT010UL-LRU	BLAUNUFT0300ULU	BLAUNUFT0200ULU	BLAUNUFT010ULU	BLAUNUFT01250UL-LRU
BLAUNUFT01250ULU	BLAUNUFT01000ULU	BLAUNUFT01000UL-LRU	BLAPM7058-AU	BLAPM7058-FU	BLAPM7059-BU
BLAPM7059-DU	BLAPM7059-FU	BLAPM7059-HU	BLAPM7059-IU	BLAPM7059-JU	BLAPM7058-JU
BLAPM7059-AU	BLAPM7058-IU	BLAPM7058-HU	BLAPS7079-HU	BLAPS7079-IU	BLAPS7079-JU
BLAPS7079-BU	BLAPS7078-JU	BLAPS7079-AU	BLAPS7078-AU	BLAPRV-10000U	BLAPRV-5000U
BLAPS7078-IU	BLAPS7078-HU	BLAPS7078-FU	BLAPS7078-DU	BLAPMP005U	BLAUNGD0201-1U
BLAUNGD0501-1U	BLAUNGD010-1U	BLAUNGD1001-1U	BLAPT7100-DU	BLAPT7101-CU	BLAPT7101-BU
BLAPT7101-DU	BLAPT7100-BU				

**Model**  
RDT3720

**UPC**  
00850084643347

**Category**  
Mixers - Fixed-Speed

**Status**  
Active