

200 series

Joystick, Z-axis control
and Touchscreen



CS200XY joystick

CS200Z Z-axis control

TS200 touchscreen with control software for Prior Scientific products

Smooth ergonomic control for stage and focus

Adjustable speed control for different applications

Operating microscope stages and focusing devices for extended periods can be strenuous. Prior Scientific's CS200 series controls are ergonomically designed according to customer feedback to ensure comfortable operation. Depending on your application, the Prior Scientific CS200XY joystick and CS200Z Z-axis controller can be used together or individually.



CS200Z Z-axis control

CS200XY joystick XY-axis control

Easy to use

The CS200Z and CS200XY have a speed control button located on the side of the unit. Separate units enable control focussing and stage simultaneously, making the microscope work faster with less user fatigue.

The **CS200XY** joystick is compatible with Prior Scientific's full range of stages.

[Motorized Stages >>](#)

Controllers:
[OptiScan ES11 >>](#)
[ProScan III Controller >>](#)

The **CS200Z** Z-axis control is compatible with Prior Scientific's full range of focus devices.

Focus products:
[Stepper Motor Focus Drive >>](#)
[Z-Axis Motorized Focusing Mounts >>](#)

Controllers:
[Controllers and Focus Controller >>](#)




TS200 touchscreen with control software for Prior Scientific products

Full system control at your fingertips


The easy-to-use TS200 touchscreen and control software is built on a large Windows touchscreen tablet. The intuitive software allows full control of all Prior Scientific products and can be used to configure product parameters to a specific application

Dimensions

CS200XY joystick



CS200Z Z-axis control



Order information

Part	Description
CS200XY	Joystick control unit
CS200Z	Focus control unit
TS200	Touchscreen control centre

Contact sales for price information >>

PRIOR[®]
S c i e n t i f i c

 prior.com
 linkedin.com/company/prior-scientific
 [@PriorScientific](https://twitter.com/PriorScientific)

UNITED KINGDOM

Prior Scientific Instruments Ltd.
Units 3-4 Fielding Industrial Estate
Wilbraham Road, Fulbourn
Cambridge, CB21 5ET
United Kingdom
Email: inquiries@prior.com
Phone: +44 (0)1223 881711

U.S.A.

Prior Scientific, Inc.
80 Reservoir Park Drive
Rockland, MA. 02370
U.S.A.
Email: info@prior.com
Phone: 781.878.8442

GERMANY

Prior Scientific Instruments GmbH
Maria Pawlowna Strasse 4
D-07743 Jena, Germany
Email: jena@prior.com
Phone: +49 (0)3641 24 20 10

JAPAN

Prior Scientific KK
Kayabacho 3rd Nagaoka Bldg 10F,
2-7-10, Nihonbashi Kayabacho, Chuo-Ku,
Tokyo103-0025, Japan
Email: info-japan@prior.com
Phone: 03-5652-8831



MicroscopeWorld

info@microscopeworld.com | 800-942-0528

CHINA

Prior Scientific Instruments (Suzhou) Ltd.
509 Tayun Plaza, 188 Tayun Road
Suzhou, 215000, China
Email: info-china@prior.com
Phone: +86 (0)512 6617 5866



Copyright © by Prior Scientific Limited, United Kingdom,
2022. Subject to changes without notice.
Version CS200 Joysticks - V2 - 1122 - EN