



DIGITAL PATHOLOGY



MoticEasyScan Two

Discover the MoticEasyScan Two, the dual-slide scanner designed for efficiency, versatility, and real-time remote evaluation.

MoticEasyScan Two is your entry point to whole slide imaging, offering real-time remote evaluation with Live View Mode, ideal for Frozen Sections and a wide range of applications. Its compact desktop design with a 2-slide capacity, optimizes lab space, while seamless multi-site connectivity facilitates effortless collaboration across labs, educational institutions, and research facilities.



Space-Efficient Design

Compact and desktop-friendly, maximizing lab space without sacrificing performance. This desktop scanner is the ideal introduction to whole-slide imaging, delivering powerful capabilities in a space-saving footprint.



2x3" Slide Compatibility

Supports optional 2x3 inch slide trays, enabling scanning of larger tissue samples for enhanced imaging versatility, accommodating diverse research and clinical needs.



Live Mode for Remote Evaluation

Operate the scanner remotely with real-time slide evaluation using Frozen Section-Live, to enable seamless collaboration. A valuable tool for medical, research, and education settings, that allows experts to connect and share insights instantly.



Flexible File Formats:

Compatible with common files such as SVS and DICOM facilitating seamless integration with third party LIS or image analysis tools.

Exceptional Image Quality:

Apochromatic objective (up to 40x, 0.13um/px) that delivers color fidelity and high-resolution whole-slide imagery.

Space-Saving Design:

Compact desktop-sized scanner tailored to fit any lab space, combining efficiency with convenience.

Capture Resolution	12 Megapixel
Resolution	0.52 $\mu\text{m}/\text{pixel}$ at 20x
	0.26 $\mu\text{m}/\text{pixel}$ at 40x
	0.13 $\mu\text{m}/\text{pixel}$ at 80x
Scanning Time	H/E stained <1 min at 20x
Objectives	10x Plan S-Apochromatic 0.3 NA
	20X Plan Apochromatic 0.75 NA
	40x Plan S-Apochromatic 0.75 NA (optional)
Slide Capacity	2
Slide Tray	2 x 2-Slide Metal Trays
Focus Mechanism	Real-time autofocus
Illumination	LED 10W
Barcode Support	1D: Interleaved 2 of 5, Code39
Length x Width x Height	205 x 392 x 420mm
Weight	12.6 kg

Contact Us:



moticrodigitalpathology.com
mdpsales@moticusa.com



info@microscopeworld.com | 800-942-0528