globe® SCIENTIFIC INC.

GCM Series

Mini Centrifuges



GCM Series Mini Centrifuges

Globe's GCM Series was designed with your lab in mind . . . whether you spin-down microcentrifuge tubes, PCR tubes or PCR strip tubes, the GCM will provide you with reliable performance use after use. Operating the GCM is quick and easy: close the lid to accelerate to 7,000rpm and simply open the lid for quick deceleration to stop.

Included with all GCM Mini Centrifuges are an eight-position microcentrifuge tube rotor for 1.5mL and 2.0mL tubes, a PCR tube rotor that accommodates 0.2mL tubes and strips as well as eight adapters for 0.2mL tubes and eight adapters for 0.5mL tubes. Changing rotors requires no tools—simply pull up to remove and press down to replace.

Features:

- Ideal for quick spin-downs and microfiltration
- EZ-Swap rotor attachments can be interchanged without tools
- Includes 8 position 1.5mL/2mL rotor, 16 position 0.2mL PCR tubes and strips rotor, 8 each 0.2ml adapters and 8 each 0.5ml adapters
- Available in your choice of 3 lid colors: red, blue and green
- Fast acceleration and braking
- "Close to start – open to stop" functionality
- 7000rpm speed (2680 x g)
- Small, lightweight and portable
- Runs quietly

Mini Centrifuges



Included accessories for GCM Series Mini Centrifuges



GCM Series Specifications	
Maximum Speed	7000rpm (fixed)
Maximum RCF	2680×g
Rotor Capacity	Eight 0.2/0.5/1.5/2.0mL tubes, sixteen 0.2mL PCR tubes or two 0.2mL PCR 8 strips
Run Time	Continuous operation
Power	20W
Noise level	≤45dB
Dimension (WxDxH)	160 x 170 x 122mm (6.3 x 6.7 x 4.8in)
Weight	0.5kg (1.1lbs)
Certifications	CE, cTÜVus, FCC, MCA
Warranty	2 Years

ORDERING INFORMATION

Item #	Description	Unit
GCM-B	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Blue Lid 120v, 60Hz, US Plug	Each
GCM-G	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Green Lid 120v, 60Hz, US Plug	Each
GCM-R	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Red Lid 120v, 60Hz, US Plug	Each
GCM-B-UK	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Blue Lid 240v, 50Hz, UK Plug	Each
GCM-G-UK	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Green Lid 240v, 50Hz, UK Plug	Each
GCM-R-UK	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Red Lid 240v, 50Hz, UK Plug	Each
GCM-B-EU	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Blue Lid 240v, 50Hz, EU Plug	Each
GCM-G-EU	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Green Lid 240v, 50Hz, EU Plug	Each
GCM-R-EU	Mini Centrifuge, 8-Place, 7000rpm Fixed Speed, Red Lid 240v, 50Hz, EU Plug	Each
GCM-ROT	Rotor, 7000rpm, 8-Place for 1.5 and 2.0mL Microcentrifuge Tubes	Each
GCM-PCR	Rotor, 7000rpm, 2 x 8 Position for 0.2mL PCR Tubes and Strips	Each
GCM-AD2	Rotor Adapter for 0.2mL Microcentrifuge Tubes, 8 each	8
GCM-AD5	Rotor Adapter for 0.5mL Microcentrifuge Tubes, 8 each	8

www.globescientific.com

Customer Service, Sales and Technical Support:
Phone 1 (800) 394-4562 • 1 (201) 599-1400
Fax 1 (201) 599-1406
E-mail mail@globescientific.com

globe® SCIENTIFIC INC.

© 2019 Globe Scientific Inc. All rights reserved.
The Globe Scientific logo is a registered trademark of Globe Scientific Inc. GCM Flyr 19064 1M